



NAD/NADH Microplate Assay Kit

User Manual

Catalog # FTA0008

(Version 2.1D)

Detection and Quantification of NAD/NADH Content in Urine,
Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and
Other biological fluids Samples.

For research use only. Not for diagnostic or therapeutic procedures.

I. INTRODUCTION.....	2
II. KIT COMPONENTS.....	3
III. MATERIALS REQUIRED BUT NOT PROVIDED.....	4
IV. SAMPLE PREPARATION.....	5
V. ASSAY PROCEDURE.....	7
VI. CALCULATION.....	8
VII. TYPICAL DATA.....	10
VIII. TECHNICAL SUPPORT.....	10
IX. NOTES.....	10

I. INTRODUCTION

Nicotinamide adenine dinucleotide (NAD^+) is a vital coenzyme found in all cells. As NAD^+ is involved in redox reactions, it is found in two forms in cells. NAD^+ is an oxidizing agent and becomes reduced to form NADH, which can be used as a reducing agent. As a result, it plays a key role in metabolism and other cellular processes. In organisms, NAD^+ can be synthesized de novo from tryptophan or aspartic acid. Because of the wide variety of functions that NAD^+ plays, it is a popular target for pharmaceuticals.

NAD/NADH Microplate Assay Kit provides a simple and direct procedure for measuring NAD^+ /NADH levels in a variety of samples. The kit is based on an alcohol dehydrogenase cycling reaction, in which the formed NADH reduces a formazan reagent. The intensity of the reduced product color, measured at 492 nm, is proportionate to the NAD^+ /NADH concentration in the sample.

II. KIT COMPONENTS

Component	Volume	Storage
96-Well Microplate	1 plate	
Assay Buffer I	30 ml x 2	4 °C
Assay Buffer II	30 ml x 2	4 °C
Reaction Buffer	10 ml x 1	4 °C
Enzyme	Powder x 1	-20 °C
Substrate	10 ml x 1	4 °C
Dye Reagent A	Powder x 1	4 °C
Dye Reagent B	1 ml x 1	4 °C
NADH Standard	Powder x 1	-20 °C, keep in dark
NAD ⁺ Standard	Powder x 1	-20 °C, keep in dark
Technical Manual	1 Manual	

Note:

Dye Reagent A: add 1 ml distilled water to dissolve before use, mix, store at 4°C.

Enzyme: add 1 ml Reaction Buffer to dissolve before use, mix, store at 4°C.

NADH Standard: add 1 ml distilled water to dissolve, mix; then add 25 µl solution into 975 µl distilled water, mix. The concentration will be 50 µmol/L.

NAD⁺ Standard: add 1 ml distilled water to dissolve, mix; then add 25 µl solution into 975 µl distilled water, mix. The concentration will be 50 µmol/L.

III. MATERIALS REQUIRED BUT NOT PROVIDED

1. Microplate reader to read absorbance at 492 nm
2. Distilled water
3. Pipettor, multi-channel pipettor
4. Pipette tips
5. Mortar
6. Ice
7. Centrifuge
8. Timer

IV. SAMPLE PREPARATION

1. For serum or plasma samples

Extract the NAD⁺:

Add 0.5 ml Assay buffer I to 0.05 ml serum or plasma; mix; boiling for 5 minutes; centrifuged at 10000g 4 °C for 10 minutes, take the supernatant into a new centrifuge tube and add 0.5 ml Assay buffer II, mix; keep it on ice for detection.

Extract the NADH:

Add 0.5 ml Assay buffer II to 0.05 ml serum or plasma; mix; boiling for 5 minutes; centrifuged at 10000g 4 °C for 10 minutes, take the supernatant into a new centrifuge tube and add 0.5 ml Assay buffer I, mix; keep it on ice for detection.

2. For tissue samples

Extract the NAD⁺:

Weigh out 0.05g tissue, homogenize with 0.5 ml Assay buffer I on ice; boiling for 5 minutes; centrifuged at 10000g 4 °C for 10 minutes, take the supernatant into a new centrifuge tube and add 0.5 ml Assay buffer II, mix; keep it on ice for detection.

Extract the NADH:

Weigh out 0.05g tissue, homogenize with 0.5 ml Assay buffer II on ice; boiling for 5 minutes; centrifuged at 10000g 4 °C for 10 minutes, take the supernatant into a new centrifuge tube and add 0.5 ml Assay buffer I, mix; keep it on ice for detection.

3. For cell and bacteria samples

Extract the NAD⁺:

Collect cell or bacteria into centrifuge tube, discard the supernatant after centrifugation, add 0.5 ml Assay buffer I for 250×10^4 cell or bacteria, sonicate (with power 20%, sonication 2s, interval 1s, repeat 30 times); boiling for 5 minutes; centrifuged at 10000g 4 °C for 10 minutes, take the supernatant into a new centrifuge tube and add 0.5 ml Assay buffer II, mix; keep it on ice for detection.

Extract the NADH:

Collect cell or bacteria into centrifuge tube, discard the supernatant after centrifugation, add 0.5 ml Assay buffer II for 250×10^4 cell or bacteria, sonicate (with power 20%, sonication 2s, interval 1s, repeat 30 times); boiling for 5 minutes; centrifuged at 10000g 4 °C for 10 minutes, take the supernatant into a new centrifuge tube and add 0.5 ml Assay buffer I, mix; keep it on ice for detection.

V. ASSAY PROCEDURE

Add following reagents into the microplate:

Reagent	Sample	Standard	Blank
Sample	20 μ l	--	--
Standard	--	20 μ l	--
Distilled water	--	--	20 μ l
Reaction Buffer	70 μ l	70 μ l	70 μ l
Enzyme	10 μ l	10 μ l	10 μ l
Substrate	80 μ l	80 μ l	80 μ l
Dye Reagent A	10 μ l	10 μ l	10 μ l
Dye Reagent B	10 μ l	10 μ l	10 μ l
Mix, keep in dark for 10 minutes at room temperature, record absorbance measured at 492 nm.			

Note:

- 1) Perform 2-fold serial dilutions of the top standards to make the standard curve.
- 2) The concentrations can vary over a wide range depending on the different samples.

For unknown samples, we recommend doing a pilot experiment & testing several doses to ensure the readings are within the standard curve range.

VI. CALCULATION

Calculation of NAD⁺:

1. According to the volume of sample

$$\begin{aligned} \text{NAD}^+ (\mu\text{mol/ml}) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / \\ &\quad V_{\text{Sample}} \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) \end{aligned}$$

2. According to the protein concentration of sample

$$\begin{aligned} \text{NAD}^+ (\mu\text{mol/mg}) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / \\ &\quad (V_{\text{Sample}} \times C_{\text{Protein}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / C_{\text{Protein}} \end{aligned}$$

3. According to the weight of sample

$$\begin{aligned} \text{NAD}^+ (\mu\text{mol/g}) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / (W \times \\ &\quad V_{\text{Sample}} / V_{\text{Assay}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / W \end{aligned}$$

4. According to the quantity of cells or bacteria

$$\begin{aligned} \text{NAD}^+ (\mu\text{mol}/10^4) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / (N \\ &\quad \times V_{\text{Sample}} / V_{\text{Assay}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / N \end{aligned}$$

Calculation of NADH:

1. According to the volume of sample

$$\begin{aligned} \text{NADH} (\mu\text{mol/ml}) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / \\ &\quad (V_{\text{Sample}} / V_{\text{Assay}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) \end{aligned}$$

2. According to the protein concentration of sample

$$\begin{aligned} \text{NADH} (\mu\text{mol/mg}) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / \\ &\quad (V_{\text{Sample}} \times C_{\text{Protein}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / C_{\text{Protein}} \end{aligned}$$

3. According to the weight of sample

$$\begin{aligned}\text{NADH } (\mu\text{mol/g}) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / (W \times \\ &\quad V_{\text{Sample}} / V_{\text{Assay}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / W\end{aligned}$$

4. According to the quantity of cells or bacteria

$$\begin{aligned}\text{NADH } (\mu\text{mol}/10^4) &= (C_{\text{Standard}} \times V_{\text{Standard}}) \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / (N \\ &\quad \times V_{\text{Sample}}/V_{\text{Assay}}) \\ &= 0.05 \times (\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}) / (\text{OD}_{\text{Standard}} - \text{OD}_{\text{Blank}}) / N\end{aligned}$$

C_{Protein} : the protein concentration, mg/ml;

C_{Standard} : the protein concentration, 50 $\mu\text{mol/L}$ = 0.05 $\mu\text{mol/ml}$;

W : the weight of sample, g;

V_{Sample} : the volume of sample, 0.02 ml;


V_{Standard} : the volume of sample, 0.02 ml;

V_{Assay} : the volume of Assay buffer, 1 ml;

N : the quantity of cell or bacteria, $N \times 10^4$.

VII. TYPICAL DATA

The standard curve is for demonstration only. A standard curve must be run with each assay.



Detection Range: 0.5 µmol/L - 50 µmol/L

VIII. TECHNICAL SUPPORT

For troubleshooting, information or assistance, please go online to www.cohesionbio.com or contact us at techsupport@cohesionbio.com

IX. NOTES