

Control Rabbit IgG Kit for Multiplex Western Blot

Catalog # ASR1003

Product Information

Description Control Rabbit IgG Kit for Multiplex Western Blot

Conjugate Unconjugated

Physical State Liquid (sterile filtered)

Buffer 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Stabilizer 50% (v/v) Glycerol

Preservative 0.01% (w/v) Sodium Azide and 0.01% (w/v) Gentamicin Sulfate

Additional Information

Shipping Condition Wet Ice

Application Note Fluorescence technology is widely used to detect proteins. However, many

common visible fluorophores often result in considerable background fluorescence in the visible range. Visible fluorophores are rarely used for membrane-based protein detection because of this high background. IRDye® 800 and IRDye® 700DX antibody and reagent conjugates are specifically designed for protein detection methods that use longer-wavelength, near-infrared (IR) fluorophores to visualize proteins in western blotting and other applications. Very low background fluorescence in the IR range provides

for a much higher signal-to-noise ratio than visible fluorophores. Detection levels in the picogram range on Western blots rival the sensitivity of chemiluminescence on film. IRDye® 800 conjugates are optimized for the Odyssey® Infrared Imaging System developed by LI-COR. IRDye® 800 conjugates are also suitable for immunofluorescence microscopy using commercially available excitation/emission filters in the 780nm/820nm range. Dual simultaneous labeling in western blots or microscopy is achieved when IRDye® 800 conjugates are used in conjunction with IRDye® 700DX or

DyLight ™680 conjugates. IRDye® 800 and IRDye® 700DX conjugates provide an ultra-sensitive and convenient alternative to standard chemiluminescent protein detection methods, as well as a valuable tool for multicolor imaging.

Storage Condition See kit insert for complete instructions.

Precautions Note This product is for research use only and is not intended for therapeutic or

diagnostic applications.

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.