

# Bovine IgG F(ab')<sub>2</sub> Peroxidase

Catalog # ASR2830

## Product Information

---

|                              |  |
|------------------------------|--|
| <b>Description</b>           | BOVINE IgG F(ab') <sub>2</sub> fragment Peroxidase conjugated          |
| <b>Conjugate</b>             | Peroxidase (Horseradish)   |
| <b>Physical State</b>        | Lyophilized  |
| <b>Host Isotype</b>          | IgG F(ab') <sub>2</sub>  |
| <b>Buffer</b>                | 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2                |
| <b>Species of Origin</b>     | Bovine   |
| <b>Reconstitution Volume</b> | 1.0 mL   |
| <b>Reconstitution Buffer</b> | Restore with deionized water (or equivalent)                           |
| <b>Stabilizer</b>            | 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free |
| <b>Preservative</b>          | None   |

## Additional Information

---

|                           |  |
|---------------------------|--|
| <b>Shipping Condition</b> | Ambient  |
| <b>Purity</b>             | This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by pepsin digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti- Bovine IgG, anti-Bovine IgG F(ab') <sub>2</sub> and anti-Bovine Serum. No reaction was observed against anti-Bovine IgG F(c) or anti-Pepsin. |
| <b>Storage Condition</b>  | Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.  |
| <b>Precautions Note</b>   | This product is for research use only and is not intended for therapeutic or diagnostic applications.  |

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.