

Anti-Cat IgM mu chain (Peroxidase Conjugated) Secondary Antibody

Rabbit Polyclonal, Peroxidase (Horseradish) Catalog # ASR2908

Product Information

Description Anti-CAT IgM mu chain (RABBIT) Antibody Peroxidase Conjugated

Host Rabbit

Conjugate Peroxidase (Horseradish)

Target Species Cat

Clonality Polyclonal Physical State Lyophilized

Host Isotype IgG

Target Isotype IgM □chain

Buffer 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Immunogen Cat IgM mu heavy chain

Reconstitution Volume 2.0 mL

Reconstitution Buffer Restore with deionized water (or equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

Preservative 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!

Additional Information

Shipping Condition Ambient

Application Note This product has been assayed against 1.0 ug of Cat IgM in a standard

capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:1,000 to 1:3,000 of the reconstitution concentration is

suggested for this product.

Purity This product is an IgG fraction antibody purified from monospecific

antiserum by a multi-step process which includes delipidation, salt

fractionation and ion exchange chromatography followed by cross adsorption

and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against

anti-Peroxidase, anti-Rabbit Serum, Cat IgM and Cat Serum. No reaction was

observed against Cat IgG.

Storage Condition Store vial at 4° C prior to restoration. For extended storage aliquot

contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted

liquid. Dilute only prior to immediate use.

Precautions NoteThis product is for research use only and is not intended for therapeutic or

diagnostic applications.

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.