

# Anti-Murine Eotaxin-2 Antibody

Catalog # ABG10088

## Product Information

---

|                    |            |
|--------------------|------------|
| <b>Application</b> | WB, E      |
| <b>Reactivity</b>  | Mouse      |
| <b>Host</b>        | Goat       |
| <b>Clonality</b>   | Polyclonal |

## Additional Information

---

|                       |                                                                                                                                                                                                                                                                                                                                                                                     |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Preparation</b>    | Produced from sera of goats pre-immunized with highly pure (>98%) recombinant mEotaxin-2. Anti-Murine Eotaxin-2 specific antibody was purified by affinity chromatography employing immobilized mEotaxin-2 matrix.                                                                                                                                                                  |
| <b>WesternBlot</b>    | To detect mEotaxin-2 by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 $\mu$ g/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant mEotaxin-2 is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.                                                                                       |
| <b>Sandwich</b>       | To detect mEotaxin-2 by sandwich ELISA (using 100 $\mu$ l/well antibody solution) a concentration of 0.5 - 2.0 $\mu$ g/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with BioGems's Biotinylated Anti-Murine Eotaxin-2 (61-065BT) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant mEotaxin-2. |
| <b>Neutralization</b> | To yield one-half maximal inhibition [ND <sub>50</sub> ] of the biological activity of mEotaxin-2 (100 ng/ml), a concentration of 7.0-8.0 $\mu$ g/ml of this antibody is required.                                                                                                                                                                                                  |
| <b>Formulation</b>    | A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution</b> | Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.                                                                                                                                                                                                                                                                                |
| <b>Storage</b>        | -20°C                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Precautions</b>    | Anti-Murine Eotaxin-2 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.                                                                                                                                                                                                                                                                    |

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.