

Cytokeratin 8 (KRT8) Antibody - With BSA and Azide

Mouse Monoclonal Antibody [Clone C-43] Catalog # AH11671

Product Information

Application IHC, IF, FC Primary Accession P05787

Other Accession 3856, 533782, 708445

Reactivity Human, Rabbit, Pig, Bovine, Sheep

Host Mouse **Clonality** Monoclonal

Isotype Mouse / IgG1, kappa

Clone Names C-43 Calculated MW 53704

Additional Information

Gene ID 3856

Other Names Keratin, type II cytoskeletal 8, Cytokeratin-8, CK-8, Keratin-8, K8, Type-II

keratin Kb8, KRT8, CYK8

Application Note IHC~~1:100~500 IF~~1:50~200 FC~~1:10~50

Storage Store at 2 to 8°C.Antibody is stable for 24 months.

Precautions Cytokeratin 8 (KRT8) Antibody - With BSA and Azide is for research use only

and not for use in diagnostic or therapeutic procedures.

Protein Information

Name KRT8

Synonyms CYK8

Function Together with KRT19, helps to link the contractile apparatus to dystrophin at

the costameres of striated muscle.

Cytoplasm. Nucleus, nucleoplasm {ECO:0000250 | UniProtKB:Q10758}.

Nucleus matrix {ECO:0000250 | UniProtKB:Q10758}

Tissue Location Observed in muscle fibers accumulating in the costameres of myoplasm at

the sarcolemma membrane in structures that contain dystrophin and spectrin. Expressed in gingival mucosa and hard palate of the oral cavity.


Background

Cytokeratin 8 (CK8) belongs to the type II (or B or basic) subfamily of high molecular weight cytokeratins and exists in combination with cytokeratin 18 (CK18). CK8 is primarily found in the non-squamous epithelia and is present in majority of adenocarcinomas and ductal carcinomas. It is absent in squamous cell carcinomas. Hepatocellular carcinomas are defined by the use of antibodies that recognize only cytokeratin 8 and 18. CK8 exists on several types of normal and neoplastic epithelia, including many ductal and glandular epithelia such as colon, stomach, small intestine, trachea, and esophagus as well as in transitional epithelium. Anti-CK8 does not react with skeletal muscle or nerve cells. Epithelioid sarcoma, chordoma, and adamantinoma show strong positivity corresponding to that of simple epithelia (with antibodies against CK8, CK18 and CK19). Reportedly, anti-CK8 is useful for the differentiation of lobular ([Pring-like, perinuclear [I]) from ductal ([Peripheral-predominant [I]) carcinoma of the breast.

References

Vojtsek et al. Folia Biol.35: 373-382 (1989). | Kovarik et al. J. Tumor Marker Oncol. 5: 219 (1990). | Bartek et al. J. Pathology,164:215-24 (1991)

Images



Formalin-fixed, paraffin-embedded human Colon carcinoma stained with Cytokeratin 8 Monoclonal Antibody (C-43).

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.