

Anti-ENTPD2 Reference Antibody (Novartis patent anti-ENTPD2)

Catalog # APR11328

Product Information


Application	Functional assay, Kinetics (SPR), Kinetics (BLI), FACS, E, FTA
Primary Accession	Q9Y5L3
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG1-REM
Calculated MW	53665

Additional Information


Target/Specificity	ENTPD2
Expression system	CHO

Protein Information


Name	ENTPD2 (HGNC:3364)
Function	Catalyzes the hydrolysis of nucleoside triphosphates (NTPs) and diphosphates (NDPs), with a marked preference for triphosphonucleosides over diphosphonucleosides (PubMed:10510450, PubMed:12888562, PubMed:16101300, PubMed:18404504). The enzyme sequentially removes phosphate groups in two successive steps, converting NTPs to nucleoside monophosphates (NMPs) via NDP intermediates (PubMed:10510450, PubMed:12888562, PubMed:16101300, PubMed:18404504). This activity contributes to the regulation of extracellular nucleotides levels (PubMed:10510450, PubMed:12888562, PubMed:16101300, PubMed:18404504). ATP hydrolysis is characterized by fast ADP accumulation and delayed AMP formation, reflecting limited ADP hydrolase activity (By similarity). Hydrolyzes ADP and UDP only to a marginal extent, and does not hydrolyze AMP (PubMed:10510450, PubMed:12888562, PubMed:16101300, PubMed:18404504). All nucleoside 5'- diphosphates are hydrolyzed at rates lower than those of their corresponding triphosphates (By similarity).
Cellular Location	[Isoform Long]: Cell membrane; Multi-pass membrane protein [Isoform gamma]: Endoplasmic reticulum membrane; Multi-pass membrane protein
Tissue Location	Brain, placenta, skeletal muscle, kidney, pancreas, heart, ovary, testis, colon, small intestine, prostate and pancreas (PubMed:9271669). No expression in adult thymus, spleen, lung, liver and peripheral blood leukocytes (PubMed:9271669)



Anti-ENTPD2 Reference Antibody (Novartis patent anti-ENTPD2) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



The purity of Anti-ENTPD2 Reference Antibody (Novartis patent anti-ENTPD2) is 97.11%, determined by SEC-HPLC.



Immobilized huNTPDase/ENTPD2-his at 8 µg/mL can bind Novartis patent anti-ENTPD2, EC₅₀= 0.01037 µg/mL.

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.