10320 Camino Santa Fe, Suite G
San Diego, CA 92121
Tel: 858.875.1900 Fax: 858.875.1999

Control Rabbit IgG Kit for Multiplex Western Blot

Catalog # ASR1003

Product Information

Description	Control Rabbit IgG Kit for Multiplex Western Blot
Conjugate	Unconjugated
Physical State	Liquid (sterile filtered)
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Stabilizer	50% (v/v) Glycerol
Preservative	0.01% (w/v) Sodium Azide and 0.01% (w/v) Gentamicin Sulfate

Additional Information

Shipping Condition	Wet Ice
Application Note	Fluorescence technology is widely used to detect proteins. However, many common visible fluorophores often result in considerable background fluorescence in the visible range. Visible fluorophores are rarely used for membrane-based protein detection because of this high background. IRDye® 800 and IRDye® 700DX antibody and reagent conjugates are specifically designed for protein detection methods that use longer-wavelength, near-infrared (IR) fluorophores to visualize proteins in western blotting and other applications. Very low background fluorescence in the IR range provides for a much higher signal-to-noise ratio than visible fluorophores. Detection levels in the picogram range on Western blots rival the sensitivity of chemiluminescence on film. IRDye® 800 conjugates are optimized for the Odyssey® Infrared Imaging System developed by LI-COR. IRDye® 800 conjugates are also suitable for immunofluorescence microscopy using commercially available excitation/emission filters in the 780nm/820nm range. Dual simultaneous labeling in western blots or microscopy is achieved when IRDye® 800 conjugates are used in conjunction with IRDye® 700DX or DyLight™680 conjugates. IRDye® 800 and IRDye® 700DX conjugates provide an ultra-sensitive and convenient alternative to standard chemiluminescent protein detection methods, as well as a valuable tool for multicolor imaging.
Storage Condition	See kit insert for complete instructions.
Precautions Note	This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.